

Product Manual

ValidorAI Hair Scanner. Clinical-grade scalp diagnostics.
Built for professionals.

Your AI hair care assistant





Product Key Features

✓ 100x magnification

✓ AI analysis backed by clinical experts


✓ 3 different light spectrums

✓ 16 paramaeters

Follicle Density, Follicle Diameter, Damage State, Hair Diameter, Scurf State, and more...

Indicator lights
Green light: the device is powered on

Power / Photo
1. Long press to turn on/off
2. Short press to capture an image



Wi-Fi
Open settings on your mobile device, go under Wi-Fi and connect to ValidorAI-XXX

How To Use

#1 App installation

iOS/Android: Search "Validor app" in the App Store or Google Play.

#2 Turn on device

Long press the power button for 3 seconds until the green light flashes.

#3 Wi-Fi connection

Connect to "ValidorAI-XXX" in your phone's Wi-Fi settings (no password).

#4 Start Detection

Open the app, tap "New scan" and follow instructions based on your scan type.



B2B SALON

Scan History

Review every supported scan for this patient, newest first.

New scan

LAST SCAN

Single Point

79

/ 100

• Good

Scalp Health

54

Follicles

Medium

Thickness

High

Density

LAST SCAN

Single Point

91

/ 100

• Excellent

Scalp Health

FAQ

❓ **Cannot connect to the device Wi-Fi?**

Ensure the device is on (green light flashing). Connect to "ValidorAI-XXX" and grant app permissions.

❓ **Wi-Fi shows no internet connection?**

This is normal. Switch to home Wi-Fi or mobile data after capturing images.

❓ **No image in the app?**

- iOS: Enable local network permissions.
- Android: Select "Yes" if prompted. Disable network acceleration if needed.

❓ **Return policy?**

14-days no-questions-asked return.

THANK YOU!

care@validorai.com • www.ValidorAI.com